

MedLED™ VISION



RIVERPOINT
MEDICAL



ZERO™ GLASSES

FOR SAFETY, VISION & LIGHT MOUNTING



Intrepid Intruder PMxtremes V2 Readers



Black Sport Frames



Multi-Color Frames



Black Frames with Cord



Built-In Magnification

MedLED Zero™ Galilean Flip-Up.
Lightweight Style.

ORDER CODES

Intrepid: Black sport frames, clear lenses _____	MLH01-Intrepid
Intruder: Colored frame; XX = Orange (OR), Black (BK), Purple (PU), Green (GR), Red (RE), Blue (BL) _____	MLH01-IntruderXX
PMxtremes: Black frames, clear lens with cord _____	MLH01-PMxtremes
V2Readers: Black frame with built-in magnification; Diopter Strength = XX = 10, 15, 20, 25, 30. _____	MLH01-V2readersXX
Cord: One spare cord (lanyard) for safety glasses _____	MLH01-Cord
Lens Cloths: Lens cleaning cloth wipes _____	MLH01-LCC
TYR Bag: Microfiber safety glasses holding bag _____	MLH01-TYRbags



ZERO™ LOUPES

BINOCULAR LOUPES



Adjustable

Knobs to adjust the intra-pupillary positioning

Optical Quality

Achromatic lenses provide a clear, crisp, and quality image. Special achromatic coating technologies reduces light reflection to 1%

MedLED Zero™ Galilean Flip-Up. Lightweight Style.

ORDER CODES



Item Code	Mag.	Working Distance	Field of View	Depth of Field
XS25340	2.5x	13-3/8" (340mm)	3-7/8" (100mm)	3-1/2" (89mm)
XS25420	2.5x	16-5/8" (420mm)	4-3/4" (120mm)	5" (127mm)
XS25460	2.5x	18-1/8" (460mm)	5-1/8" (130mm)	3-1/2" (89mm)
XS25500	2.5x	19-5/8" (500mm)	5-1/2" (140mm)	5" (127mm)
XS30340	3.0x	13-3/8" (340mm)	2-1/2" (65mm)	3" (77mm)
XS30420	3.0x	16-5/8" (420mm)	3-1/2" (90mm)	6" (152mm)
XS30460	3.0x	18-1/8" (460mm)	4-3/8" (110mm)	5-7/8" (150mm)
XS30500	3.0x	19-5/8" (500mm)	5" (125mm)	6-1/2" (165mm)
XS35340	3.5x	13-3/8" (340mm)	2-3/8" (60mm)	2-1/2" (65mm)
XS35420	3.5x	16-5/8" (420mm)	3-1/2" (90mm)	4-3/8" (110mm)
XS35460	3.5x	18-1/8" (460mm)	3-7/8" (100mm)	5-1/8" (130mm)
XS35500	3.5x	19-5/8" (500mm)	4-1/2" (115mm)	6-7/8" (175mm)

MedLED Zero™ Keplerian Flip-Up. Higher Magnification for Advanced Procedures.



Item Code	Mag.	Working Distance	Field of View	Depth of Field
XK40340	4.0x	13-3/8" (340mm)	2-1/8" (55mm)	12-14" (300-365mm)
XK40420	4.0x	16-5/8" (420mm)	2-3/8" (60mm)	14-18" (360-460mm)
XK40460	4.0x	18-1/8" (460mm)	2-1/2" (62mm)	16-21" (405-530mm)
XK40500	4.0x	19-5/8" (500mm)	2-3/4" (70mm)	17-23" (425-580mm)
XK45340	4.5x	13-3/8" (340mm)	1-3/4" (45mm)	12-15" (310-370mm)
XK45420	4.5x	16-5/8" (420mm)	2" (50mm)	15-18" (380-465mm)
XK45460	4.5x	18-1/8" (460mm)	2-3/8" (60mm)	17-20" (440-510mm)
XK45500	4.5x	19-5/8" (500mm)	2-1/2" (65mm)	17-22" (440-565mm)

MANUFACTURED
IN PORTLAND,
OREGON USA



Med**LED**[™]

BEING THE FIRST AND STILL THE BEST IN SURGICAL LED LIGHTING

Mission. Freedom without sacrifice. Providing hi-tech lighting alternatives to surgeons looking to be free of cables and expensive maintenance costs without sacrificing basic surgical lighting needs.



MADE IN THE U.S.A.

A Brand Under Riverpoint Medical
825 NE 25th Ave., Portland, OR 97232, USA
WEB www.rpmed.com/medled
TOLL FREE +01 (866) 445-4923
ISO 13485:2003 | ISO 9001:2008 | FDA Registered

PUB159-01 10-2013