

RP ONCOLOGY

BRACHYTHERAPY NEEDLES

BRACHYTHERAPY ACCESSORIES

GOLD FIDUCIAL MARKERS



RP ONCOLOGY

BRACHYTHERAPY NEEDLES
BRACHYTHERAPY ACCESSORIES
GOLD FIDUCIAL MARKERS



BRACHYTHERAPY NEEDLES


ROUND HUB NEEDLE

RP-1120-1820.....18g (1.2mm x 20cm); Blue foil marker / Blue stilet hub



SQUARE HUB NEEDLE

RP-1100-1820.....18g (1.2mm x 20cm); Blue foil marker / Blue stilet hub
RP-1100Y-1820.....18g (1.2mm x 20cm); Blue foil marker / Yellow stilet hub
RP-1100-1812.....18g (1.2mm x 12cm); Blue foil marker / Blue stilet hub




ROUND HUB

SQUARE HUB

- Electropolished for sharpness
- Siliconized for easy insertion
- Mirror polished I.D. for a smooth, easy seed deployment
- Echo-enhanced tipped for ultrasound visualization
- Centimeter markings for depth placement
- Rounded stilet for better seed delivery
- Foil hub marker to indicate bevel tip orientation

APPLICATOR NEEDLE

RP-1200-1820.....18g (1.2mm x 20cm)
RP-1200-1720.....18g (1.4mm x 20cm)



APPLICATOR

- Electropolished for sharpness and siliconized for easy insertion
- Mirror polished I.D. for a smooth, easy seed deployment
- Echo-enhanced tipped for ultrasound visualization
- Centimeter markings for depth placement
- Trocar stilet for smooth direct entry
- Designed to fit the Mick® 200 T.P.V. Applicator

BIOPSY NEEDLES



MAGNUM STYLE




RP-1000-1820.....18g (1.2mm x 20cm); Blue foil marker / Blue stylet hub, Round Hub
 RP-1000-1825.....18g (1.2mm x 25cm); Blue foil marker / Blue stylet hub, Round Hub

BIOPTY / PRO-MAG* STYLE



RPM-1120-1820.....18g (1.2mm x 20cm); Blue foil marker / Blue stylet hub, Round Hub



- Etched Cutting notch for better tissue retention
- Echo enhanced tip for ultrasound placement
- Spacer clip for easy loading
- Centimeter markings for depth placement
- Electropolished for sharpness & smooth entry
- 19mm or 17mm cutting notch for tissue capture

*BIOPTY / PRO-MAG

MAGNUM

BRACHYTHERAPY ACCESSORIES

GRID

RPG-18....18 Guage Grid



- Easy Needle Insertion
- Precise 18g Needle Channel Guides
- Accurate Placement
- Single, One-Time Sterile Use
- Ensure that Ultrasound Tip Matches Grid Graphics



GOLD FIDUCIAL MARKERS

KITS

GM 1231-1720.....(1) 1.2 x 3mm Gold Marker w/ 17g x 20cm Needle
 GM 1233-1720.....(3) 1.2 x 3mm Gold Marker w/ 17g x 20cm Needle
 GM 1251-1720.....(1) 1.2 x 5mm Gold Marker w/ 17g x 20cm Needle
 GM 1253-1720.....(3) 1.2 x 5mm Gold Marker w/ 17g x 20cm Needle
 GM 1031-1820.....(1) 1.0 x 3mm Gold Marker w/ 18g x 20cm Needle
 GM 1033-1820.....(3) 1.0 x 3mm Gold Marker w/ 18g x 20cm Needle
 GM 1051-1820.....(1) 1.0 x 5mm Gold Marker w/ 18g x 20cm Needle
 GM 1053-1820.....(3) 1.0 x 5mm Gold Marker w/ 18g x 20cm Needle
 GM 0833-1820.....(1) 0.8 x 3mm Gold Marker w/ 18g x 20cm Needle
 GM 0851-1820.....(1) 0.8 x 5mm Gold Marker w/ 18g x 20cm Needle



*Custom needle molding and manufacturing upon request.

*BIOPTY® is a registered trademark of BARD®

*PRO-MAG™ is a trademark of ANGIOTECH



MANUFACTURED
IN PORTLAND,
OREGON USA

RP ONCOLOGY.

Riverpoint Medical designs, develops, and markets next generation unique medical devices for the veterinary and healthcare markets.



MADE IN THE U.S.A.

825 NE 25th Ave., Portland, OR 97232, USA

WEB www.rpmed.com

EMAIL info@rpmed.com

OFFICE +01 (503) 517-8001

TOLL FREE +01 (866) 445-4923

FAX +01 (503) 517-8002

ISO 13485:2003 | ISO 9001:2008 | FDA Registered