

# RiverBand™

n-Butyl Cyanoacrylate  
SURGICAL ADHESIVE

# RIVERPOINT MEDICAL

**RiverBand™** contains purified n-Butyl Cyanoacrylate with a D&C certified dye that allows for differentiation between the product and bodily fluid.

Surgical Adhesive should be applied topically; it is not intended for internal closures.

**RiverBand™** Surgical Adhesive **dries clear** for clean appearance. It **sets up quickly** (5-15 seconds) and **dries flexible** for maintaining closure on wound site.

The included disposable applicator tip fits snugly into the tip on the bottle and allows greater control over the size and placement to the drop of product.



**RiverBand™** is designed for wound management in situations such as:

- lacerations and abrasions
- deglaws
- cuts
- tail docking

Package includes:  
2.0 mL dropper bottle  
12 disposable applicator tips  
package insert.

**RIVERPOINT  
MEDICAL**

825 NE 25th Avenue  
Portland, Oregon 97232 USA  
Office: 503 517-8001  
Toll Free: 866 445-4923  
Fax: 503 517-8002  
[www.rpmed.com](http://www.rpmed.com)