


RP Cutting Biopsy Needle

**RIVERPOINT
MEDICAL**



- Etched Cutting notch for better tissue retention
- Echo enhanced tip for ultrasound placement
- Spacer clip for easy loading

- Centimeter markings for depth placement
- Electropolished for sharpness & smooth entry
- 19mm or 17mm cutting notch for tissue capture

Magnum Style

Biopty / Promag



Cutting notch

Order Codes:

Magnum Style
RP 1820 18g x 20cm
RP 1825 18g x 25cm
Biopty / Promag
RP 1820m 18g x 20cm

**RIVERPOINT
MEDICAL**

Riverpoint Medical
825 NE 25th Avenue
Portland, OR 97232
(866) 445-4923