

# RP Sleeves and Spacers

*Innovators In Brachytherapy Since 1998*

## Rapidsorb Multifil™ Carrier

RP Synthetic Absorbable PGLA 90/10 Braided Sleeve Design



RPBS-001

\*FDA Registered

## Rapidsorb Mono™ Carrier

RP Synthetic Absorbable PGLA 20/80 Hollow Mono Sleeve Design



RPMS-001

\*FDA Registered

## RP Precision™ Cut Brachy-Spacer

Synthetic Absorbable 90/10 PGLA

**RPS-910-50-50**

PGLA 90/10 5.0mm spacer qty: 50

**RPS-910-55-150**

PGLA 90/10 5.5mm spacer qty: 150

**RPS-910-100-25**

PGLA 90/10 1.0cm spacer qty: 25



**Available in 5.0mm / 5.5mm / 1cm lengths**

**Custom lengths available upon request**

\*Non-Sterile or Sterile

**RIVERPOINT  
MEDICAL**



Riverpoint Medical  
825 NE 25th Avenue  
Portland, OR 97232  
(866) 445-4923