



RIVERPOINT
MEDICAL



SUTURE BROCHURE

Absorbable Suture

Vilet™
Coated Braided Violet
Polyglactin 910
Synthetic Absorbable Suture
CE MARK, AVAILABLE FOR EXPORT ONLY
Equivalent to Vicryl*

ViletQ™
fast absorbing
Equivalent to Vicryl Rapide*



Vilet™ PLUS
Antibacterial
Coated Braided
Polyglactin 910
Synthetic Absorbable Suture
CE MARK, AVAILABLE FOR EXPORT ONLY
Equivalent to Vicryl PLUS*
VET USE ONLY




Chromic Gut
Chromic Treated
Natural Absorbable Suture
FOR VET USE ONLY



PDM™
Violet Monofilament
Polydioxanone
Synthetic Absorbable Suture
Equivalent to PDS*II



Mono Q™
Violet Monofilament
POLY(glycolide-co-caprolactone) PGCL
Synthetic Absorbable Suture
Equivalent to Monocryl*



Vilet™ II
Coated Braided Violet
100% Polyglycolic Acid
Synthetic Absorbable Suture
AVAILABLE IN U.S. ONLY



Nonabsorbable Suture


Riverlon™
Blue or Black
Nylon Monofilament
Nonabsorbable Suture
Equivalent to Ethilon*



Rivermid™
Braided Twisted Nylon
Nonabsorbable Suture
Equivalent to Supramid**



HS FIBER™
Ultra High Molecular Weight Polyethylene
Nonabsorbable Suture



Riverpro™
Blue
Polypropylene Monofilament
Nonabsorbable Suture
Equivalent to Prolene*



Fluoropro™
Pink Fluorescent
Polypropylene Monofilament
Nonabsorbable Suture
Equivalent to Fluorofil***



RP NEEDLES

3/8 Circle Reverse Cutting		3/8 Circle Hygienic Reverse Cutting	
	FS-2 19mm		PS-2 19mm
	FS-1 24mm		PS-1 24mm
	FS 26mm		
	FSL 30mm		
	FSLX 40mm		
1/2 Circle Reverse Cutting		1/2 Circle Taper Point	
	X-1 22mm		RB-1 17mm
	CP-2 26mm		SH-1 22mm
	CP-1 36mm		SH 26mm
	CP 40mm		MO-6 26mm
			CT-3 22mm
			CT-2 26mm
			CT-1 36mm
Spitula		1/2 Circle Precision Taper	
	CS160-8 5.5mm		V-34 36mm
	S-29 8.0mm		
	TG140-8 6.5mm		


Straight Reverse Cutting	
	KS 60mm

Additional needles include, LR, TP-1, CCS-1 & CCS.
(Other needles available upon request)

the 300 Difference


Ductile Strength

100% increase in ductility.
When needle bends it is less likely to break when reshaping.



Sharper Needle

12% increase in penetration.
Less force, less trauma.



*registered trademark of ETHICON, INC

**registered trademark of S. Jackson, INC

***Fluorofil is a registered trademark of Schering Plough Animal Health Corp.

RP SUTURE PRODUCTS

PDM™ (polydioxanone)

RP CODE	ETHICON CODE	NEEDLE	SIZE		LENGTH	
			U.S.P.	EP		
D304	Z304H	RB-1	4-0	1.5	30"	75cm
D305	Z305H	RB-1	3-0	2	30"	75cm
D310	Z310H	SH-1	4-0	1.5	30"	75cm
D311	Z311H	SH-1	3-0	2	30"	75cm
D315	Z315H	SH	4-0	1.5	30"	75cm
D316	Z316H	SH	3-0	2	30"	75cm
D317	Z317H	SH	2-0	3	30"	75cm
D332	Z332H	CT-2	3-0	2	30"	75cm
D333	Z333H	CT-2	2-0	3	30"	75cm
D334	Z334H	CT-2	0	3.5	30"	75cm
D338	Z338H	CT-1	3-0	2	30"	75cm
D339	Z339H	CT-1	2-0	3	30"	75cm
D340	Z340H	CT-1	0	3.5	30"	75cm
D341	Z341H	CT-1	1	4	30"	75cm
D397	Z397H	FS-2	4-0	1.5	30"	75cm
D398	Z398H	FS-2	3-0	2	30"	75cm
D451	Z451H	FS-1	2-0	3	30"	75cm
D452	Z452H	FS-1	3-0	2	30"	75cm
D466	Z466H	CP-1	2-0	3	30"	75cm
D467	Z467H	CP-1	0	3.5	30"	75cm
D468	Z468H	CP-1	1	4	30"	75cm
D968	Z968H	CP-2	3-0	2	30"	75cm
D969	Z969H	CP-2	2-0	3	30"	75cm
D970	Z970H	CP-2	0	3.5	30"	75cm

Mono Q™ (poliglecaprone 25)

RP CODE	ETHICON CODE	NEEDLE	SIZE		LENGTH	
			U.S.P.	EP		
Q315	Y315H	SH	4-0	1.5	30"	75cm
Q316	Y316H	SH	3-0	2	30"	75cm
Q317	Y317H	SH	2-0	3	30"	75cm
Q344	Y344H	CT-1	3-0	2	36"	90cm
Q345	Y345H	CT-1	2-0	3	36"	90cm
Q346	Y346H	CT-1	0	3.5	36"	90cm
Q463	Y463G	P-3	5-0	1	18"	45cm
Q464	Y464G	P-3	4-0	1.5	18"	45cm
Q732	Y732H	SH	2-0	3	36"	90cm
Q762	Y762H	CT-2	2-0	3	36"	90cm
Q922	Y922H	FS-2	4-0	1.5	36"	90cm
Q923	Y923H	FS-2	3-0	2	36"	90cm
Q942	Y942H	FS-1	3-0	2	36"	90cm
Q943	Y943H	FS-1	2-0	3	36"	90cm
Q967	Y967G	CP-1	0	3.5	36"	90cm
Q987	Y987G	FSL	0	3.5	36"	90cm

Riverlon™ (blue nylon)

RP CODE	ETHICON CODE	NEEDLE	SIZE		LENGTH	
			U.S.P.	EP		
460BL	460H	LR	2-0	3	20"	50cm
628BL	628H	KS	2-0	3	30"	75cm
662BL	662H	FS-2	4-0	1.5	30"	75cm
663BL	663H	FS-1	3-0	2	30"	75cm
664BL	664H	FS	2-0	3	30"	75cm
665BL	665H	FSL	0	3.5	30"	75cm
1674BL	1674H	FSLX	2-0	3	30"	75cm

Riverlon™ (black nylon)

RP CODE	ETHICON CODE	NEEDLE	SIZE		LENGTH	
			U.S.P.	EP		
662BK	662H	FS-2	4-0	1.5	18"	45cm
663BK	663H	FS-1	3-0	2	18"	45cm
664BK	664H	FS	2-0	3	18"	45cm

Vilet™ (polyglactin 910)

RP CODE	RP CODE	ETHICON CODE*	NEEDLE	SIZE		LENGTH	
				U.S.P.	EP		
V195	VG195	J195H	CP	2	5	30"	75cm
V311	VG311	J311H	SH	3-0	2	30"	75cm
V315	VG315	J315H	SH	4-0	1.5	30"	75cm
V316	VG316	J316H	SH	3-0	2	30"	75cm
V317	VG317	J317H	SH	2-0	3	30"	75cm
V327	VG327	J327H	CT-3	3-0	2	30"	75cm
V332	VG332	J332H	CT-2	3-0	2	30"	75cm
V333	VG333	J333H	CT-2	2-0	3	30"	75cm
V334	VG334	J334H	CT-2	0	3.5	30"	75cm
V338	VG338	J338H	CT-1	3-0	2	30"	75cm
V339	VG339	J339H	CT-1	2-0	3	30"	75cm
V340	VG340	J340H	CT-1	0	3.5	30"	75cm
V341	VG341	J341H	CT-1	1	4	30"	75cm
V345	VG345	J345H	CT-1	2-0	3	36"	90cm
V397	VG397	J397H	FS-2	4-0	1.5	30"	75cm
V398	VG398	J398H	FS-2	3-0	2	30"	75cm
V452	VG452	J452H	FS-1	3-0	2	30"	75cm
V453	VG453	J453H	FS-1	2-0	3	30"	75cm
V460	VG460	J460H	X-1	3-0	2	30"	75cm
V461	VG461	J461H	X-1	2-0	3	30"	75cm
V466	VG466	J466H	CP-1	2-0	3	30"	75cm
V467	VG467	J467H	CP-1	0	3.5	30"	75cm
V468	VG468	J468H	CP-1	1	4	30"	75cm
V586	VG586	J586H	FSL	2-0	3	30"	75cm
V587	VG587	J587H	FSL	0	3.5	30"	75cm
V808	VG808	J808H	CCS-1	3	6	30"	75cm
V880	VG880	J880H	TP-1	2	5	54"	135mm

* For ETHICON PLUS Code change prefix "J" to "VCP"

Chromic Gut

RP CODE	ETHICON CODE	NEEDLE	SIZE		LENGTH	
			U.S.P.	EP		
CG122	122H	SH	3-0	3	30"	75cm
CG123	123H	SH	2-0	3.5	30"	75cm
CG636	636H	FS-2	3-0	3	30"	75cm
CG637	637H	FS-1	2-0	3.5	30"	75cm
CG812	812H	CT-1	0	4	30"	75cm
CG882	882H	CT-2	3-0	3	30"	75cm
CG883	883H	CT-2	2-0	3.3	30"	75cm

Fluoropro™ (polypropylene)

RP CODE	NEEDLE	SIZE		LENGTH	
		U.S.P.	EP		
F8665	FS-2	3-0	2	30"	75cm
F8683	FS-2	4-0	1.5	30"	75cm
F8684	FS-1	3-0	2	30"	75cm
F8685	FS	2-0	3	18"	45cm
F8689	FSLX	2-0	3	30"	75cm
F8690	FSLX	0	3.5	30"	75cm

Riverpro™ (polypropylene)

RP CODE	ETHICON CODE	NEEDLE	SIZE		LENGTH	
			U.S.P.	EP		
P8411	8411H	CT-2	2-0	3	30"	75cm
P8424	8424H	CT-1	0	3.5	30"	75cm
P8444	8444H	V-34	0	3.5	30"	75cm
P8665	8665H	FS-2	3-0	2	18"	45cm
P8684	8684H	FS-1	3-0	2	18"	45cm
P8685	8685H	FS	2-0	3	18"	45cm
P8690	8690H	FSLX	0	3.5	30"	75cm
P8825	8825H	TP-1	2	5	60"	150cm
P8872	8872H	RB-1	3-0	2	30"	75cm

(Other codes may be available upon request)

RP SUTURE CASSETTE PRODUCTS

(VET USE ONLY)

PDM™ (polydioxanone)

RP CODE	SUTURE SIZE	LENGTH
MCD-2	2	25 m
MCD-1	1	25 m
MCD-0	0	25 m
MCD-2/0	2-0	25 m
MCD-3/0	3-0	25 m
MCD-4/0	4-0	25 m

Mono Q™ (poliglecaprone 25)

RP CODE	SUTURE SIZE	LENGTH
MCQ-2	2	25 m
MCQ-1	1	25 m
MCQ-0	0	25 m
MCQ-2/0	2-0	25 m
MCQ-3/0	3-0	25 m
MCQ-4/0	4-0	25 m


Chromic Gut

RP CODE	SUTURE SIZE	LENGTH
MCG-4	4	25 m
MCG-3	3	40 m
MCG-2	2	45 m
MCG-1	1	50 m
MCG-0	0	50 m
MCG-2/0	2-0	50 m
MCG-3/0	3-0	50 m
MCG-4/0	4-0	50 m

Fluoropro™ (polypropylene)

RP CODE	SUTURE SIZE	LENGTH
MCF-0	0	85 m
MCF-2/0	2-0	110 m
MCF-3/0	3-0	110 m
MCF-4/0	4-0	110 m

Absorption chart



Vilet™ / Vilet PLUS™ (polyglactin 910)

RP (Vilet) CODE	RP (Vilet Plus) CODE	SUTURE SIZE	LENGTH
MCV-2	MCVP-2	2	25 m
MCV-1	MCVP-1	1	25 m
MCV-0	MCVP-0	0	25 m
MCV-2/0	MCVP-2/0	2-0	25 m
MCV-3/0	MCVP-3/0	3-0	25 m
MCV-4/0	MCVP-4/0	4-0	25 m

Rivermid™ (coated twisted nylon)

RP CODE	SUTURE SIZE	LENGTH
MCS-4	4	25 m
MCS-3	3	50 m
MCS-2	2	50 m
MCS-1	1	75 m
MCS-0	0	75 m
MCS-2/0	2-0	100 m
MCS-3/0	3-0	100 m
MCS-4/0	4-0	100 m

Riverlon™ (nylon)

RP CODE	SUTURE SIZE	LENGTH
MCN-3	3	50 m
MCN-2	2	50 m
MCN-1	1	75 m
MCN-0	0	75 m
MCN-2/0	2-0	100 m
MCN-3/0	3-0	100 m



RiverBand™ Tissue Adhesive

RP CODE	SIZE	MATERIAL	APPLICATORS
TA-01	2mL	n-Butyl Cyanoacrylate	12 Applicator tips



825 NE 25th Avenue
 Portland, Oregon 97232 USA
 Office: 503 517-8001
 Toll Free: 866-445-4923
 Fax: 503 517-8002
 info@rpmed.com

www.rpmed.com

Registered trademarks of ETHICON, INC.
 *Registered trademark of S. Jackson, INC.
 **Fluorod is a registered trademark of Scoring Plough Animal Health Corp.